

A broad spectrum bacteriostatic and bactericidal activity against gram negative and gram-positive organisms

Catalog	Unit Size
TBS8037-50ML	50 mL
TBS8037-100ML	100 mL

DESCRIPTION

Penicillin-Streptomycin-Fungizone (Amphotericin B) is a broad spectrum bacteriostatic and bactericidal, with activity against gram negative and gram-positive organisms including yeast and mold.

COMPONENTS:

Fungizone (Amphotericin B): 25 µg/ml

Potassium Penicillin: 10,000 units/ml

Streptomycin Sulfate: 10,000 µg/ml

NaCl: 8.5 g/L

PACK SIZE

50mL or 100 mL /bottle or 1x1 L/bottle, Sterile.

Store at -20°C.

Research use only.

RELATED PRODUCTS

MSC Medium (TBS8021)

Ham's F-12 Medium (TBS8032)

DMEM Medium (TBS8061)

Chondrogenic Differentiation Medium (TBS8062)

RPMI-1640 Medium (TBS8063)

Hybridoma Growth Medium (TBS8074)

DMEM-F12 Mixture (TBS8083)

ESC/iPSC-qualified FBS (TBS8002)

Adipocyte Differentiation Cocktail (TBS8017)

0.1% Gelatin Solution (TBS8004)

2x HBS, pH7.05 (TBS5076)

Cell Culture Grad Water (TBS5050) LB Medium (TBS8056)

SOB Medium (TBS8057)

SOC Broth Medium (TBS8058)

2xYT Broth Medium (TBS8059)