

Anti-Canine Procalcitonin (PCT) Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TMC1023-0.2MG	0.2mg
TMC1023-1MG	1mg

Description: Anti-canine procalcitonin (PCT) Mab to canine PCT. This item (TMC1023) can be paired with TMC1024 for ELISA and Lateral Flow.

Specificity:	This antibody reacts with canine/dog PCT.
Host Animal:	Balb/c mouse.
Isotype:	IgG2b.
Source:	Mouse ascites.
Immunogen:	Canine PCT protein.
Format:	Purified, liquid.
Purification:	Protein A Purified Monoclonal Antibody.
Concentration:	1.0 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TMC1024 for canine/dog PCT detection.
Application:	For canine/dog PCT quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	Store at -20°C for 2 years.

For Research Use only