

## Benzamil HCl, Sodium channel blocker

---

Catalog	Unit
TBI4541-10MG	10 mg
TBI4541-50MG	50 mg

### Product Details

**Formal Name:** N-(Benzylamidino)-3,5-diamino-6-chloropyrazinecarboxamide hydrochloride

**Molecular Formula:** C<sub>13</sub>H<sub>14</sub>ClN<sub>7</sub>O HCl

**Formula Weight:** 356.21

**CAS Number:** 161804-20-2


**Purity:** >98%

**Formulation:** powder

**Solubility:** Soluble in DMSO (up to 25 mg/ml) or in Water (up to 5 mg/ml).

**Storage:** -20°C

**Stability:** ≥ 2 years.



### Applications

Sodium channel blocker

### Functions

Sodium channel blocker

### Application Procedures

First dissolved in DMSO (up to 25 mg/ml) or in Water (up to 5 mg/ml), then diluted to aqueous buffer. Solutions in DMSO or Water may be stored at -20°C for up to 3 months.

**For research use only.**