

PDE2A1 / Phosphodiesterase 2A1 Antibody (TBS10040)

Specification Sheet

Catalog Number	Unit Size
TBS10040-100UL	100 μ L
TBS10040-200UL	200 μ L

Description: Anti-phosphodiesterase 2A1 (PDE2A1) polyclonal antibody N-epitope. This item (TBS10040) is validated for ELISA, Western Blot, IHC, IP and other immunoassays.

Specificity: This antibody only reacts with PDE2A1

Host Animal: Rabbit.

Isotype: IgG.

Clonality: Polyclonal.

Determinant: N-epitope.

Species Cross Reactivity: Human, Monkey, Mouse, Rat.

Immunogen: Synthetic peptide corresponding to region within amino acids 1-50 on N-terminus of human PDE2A1 protein.

Format: Affinity purified, liquid.

Purification: Affinity purified.

Concentration: 0.5mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 0.05% NaN_3 , pH 7.2.

Pairing Assays: No.

Application:

Assay Method	Dilution
ELISA	1:10,000
Immunohistochemistry	1:100
Immunoprecipitation	1:200
Western Blot	1:500

Storage: -20°C for 2 years.

For Research Use only