

PDE2A3 / Phosphodiesterase 2A3 Antibody (TBS10042)

Specification Sheet

Catalog Number	Unit Size
TBS10042-100UL	100 µL
TBS10042-200UL	200 µL

Description: Anti-phosphodiesterase 2A3 (PDE2A3) polyclonal antibody N-epitope. This item (TBS10042) is validated for ELISA, Western Blot, IHC, IP and other immunoassays.

Specificity: This antibody only reacts with PDE2A3.

Host Animal: Rabbit.

Isotype: IgG.

Clonality: Polyclonal.

Determinant: N-epitope.

Species Cross Reactivity: Human, Monkey, Mouse, Rat.

Immunogen: Synthetic peptide corresponding to region within amino acids 1-50 on N-terminus of human PDE2A2 protein.

Format: Affinity purified, liquid.

Purification: Affinity purified.

Concentration: 0.5mg/ml.

Buffer: 10 mM PBS, 0.05% NaN₃, pH 7.2.

Pairing Assays: No.

Application:

Assay Method	Dilution
ELISA	1:10,000
Immunohistochemistry	1:100
Immunoprecipitation	1:200
Western Blot	1:500

Storage: -20°C for 2 years.

For Research Use only