

Optical 384-well QPCR Plate

Catalog	Unit
L0008-10	10 Plates / Box
L0008-50	50 Plates / Case

Product Details

Materials: Polypropylene.

Product Type: Microplate, DNA/RNA/RNase Free.

Number of Wells: 384.

Volume of Well: 40 μ L.

Skirt Style: Skirted.

Pack Size: 10 plates/ box.

Case Size: 50 plates (5 boxes).

Application: Use for qPCR amplification in 384 well qPCR machine System including Applied Biosystems Real-Time PCR Systems, GeneAmp PCR System, and other qPCR machine systems. The plate is compatible with MicroAmp™ Optical 384-Well Reaction Plate (4309849).

Reaction Speed: Fast or Standard.

Storage: Room Temperature.

For research use only.