

Phospho-Tau Thr181 (p-tau181) Mouse Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10034-0.1MG	0.1mg
TBS10034-0.5MG	0.5mg

Description:	Phospho-Tau Thr181 (p-tau181) Mouse Monoclonal Antibody. This antibody is validated as detection paired with TBS10193 (capture) in ELISA. It can be used for Western Blot, IHC, lateral flow, and chemiluminescent.
Specificity:	This antibody reacts with synthesized p-tau181 protein.
Host Animal:	Balb/c mouse.
Isotype:	IgG1.
Immunogen:	P-tau181 peptide.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	1.0 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH7.2.
Pairing Assays:	This antibody as detection is paired with TBS10193 as capture.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	ELISA, Western Blot, and IHC, lateral flow, and chemiluminescent.
Storage:	-20°C for 2 years.

For Research Use only