

Tau Mouse Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10049-0.1MG	0.1mg
TBS10049-0.5MG	0.5mg

Description:	Tau Mouse Monoclonal Antibody. This antibody can recognize phosphor tau and non-phosphor tau proteins. It can be used with other tau antibodies for ELISA. It can be used for Western Blot, and IHC.
Specificity:	This antibody reacts with phosphor tau or non-phosphor proteins in human, mouse, and rats.
Species Reactivity:	Hu, Ms, and Rt
Host Animal:	Balb/c mouse.
Isotype:	IgG1.
Immunogen:	Tau peptide.
Format:	Purified, liquid.
Purification:	Chromatography.
Concentration:	1.0 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH7.2.
Pairing Assays:	Not determined.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	ELISA, Western Blot, and IHC.
Storage:	-20°C for 3 years.

For Research Use only