

Tau 441 Rabbit Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10050-0.1MG	0.1mg
TBS10050-0.5MG	0.5mg

Description:	Tau 441 Rabbit Monoclonal Antibody. This antibody can recognize phosphor tau and non-phosphor tau proteins. It can be used with other tau antibodies for ELISA. This antibody can be paired with Tau Mouse monoclonal antibody (TBS10049) for Total tau ELISA detection. It can be used for Western Blot, and IHC.
Specificity:	This antibody reacts with phosphor tau or non-phosphor proteins in human, mouse, and rats.
Species Reactivity:	Hu, Ms, and Rt
Host Animal:	Rabbit.
Isotype:	IgG1.
Immunogen:	Tau peptide.
Format:	Purified, liquid.
Purification:	Affinity Chromatography.
Concentration:	0.5 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH7.2.
Pairing Assays:	Not determined.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	ELISA, Western Blot, and IHC.
Storage:	-20°C for 3 years.

For Research Use only