

F-actin Mouse Monoclonal Antibody

Catalog	Unit
TBS10086-100UG	100 ug
TBS10086-500UG	500 ug

Product Details

Description: F-actin is involved in numerous cellular processes including cell motility, cytokinesis, vesicle and organelle transport, muscle contraction, and maintaining cell morphology.

Applications: WB, IHC, ICC, ELISA

Host: Mouse

Species Cross Reactivity: Human, Mouse, Rat

Isotype: IgM


Concentration: 1 mg/mL

Formulation: PBS, pH 7.4; 0.09% sodium azide

Clonality: Monoclonal

Storage: 4°C

Images



Western Blotting (dilution range 1:100-1:500), Immunocytochemistry (dilution range 1:10), ELISA (dilution range 1:10).

For research use only.