

Phospho-Tau Thr181 (p-tau181) Rabbit Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10193-0.1MG	0.1mg
TBS10193-0.5MG	0.5mg

Description:	Phospho-Tau Thr181 (p-tau181) Rabbit Monoclonal Antibody. This antibody is validated as capture paired with TBS10034 (detection) in ELISA. It can be used for Western Blot, IHC, Lateral Flow, and chemiluminescent.
Specificity:	This antibody reacts with synthesized p-tau181 protein.
Host Animal:	Rabbit.
Isotype:	IgG1.
Species Reactivity:	Human, Mouse, Rat.
Immunogen:	P-tau181 peptide.
Format:	Purified, liquid.
Purification:	Affinity purification.
Concentration:	1.0 mg/ml.
Affinity Constant:	Not determined.
Formulation:	PBS with 0.02% sodium azide, pH7.3.
Pairing Assays:	This antibody can be paired with TBS10034 in ELISA, Lateral Flow, and chemiluminescent.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	ELISA, Western Blot, IHC and Lateral Flow, and chemiluminescent.
Storage:	-20°C for 2 years.

For Research Use only