

<b>Catalog</b>	<b>Unit</b>
TBP0169-0.1	0.1 mg
TBP0169-1	1 mg
TBP0169-5	5 mg
TBP0169-10	10 mg

**Description**

Tribioscience Recombinant Human Insulin-like Growth Factor (rhIGF-1) is manufactured through a recombinant expression system, guaranteeing consistent quality and high purity. This animal-free protein, with its low endotoxin levels, offers superior safety and reliability for research use.

**Specification**

**Source:** Escherichia coli

**Appearance:** White lyophilized powder

**Molecular weight:** 8.3 kDa

**Formulation:** It is lyophilized with mannitol as the stabilizer.

**Purity:** More than 95% as determined by SDS-PAGE.


**Endotoxin:** Less than 5 EU/mg

**Storage:** Stored at temperatures -20°C for at least 36 months.

**Applications:** Cell culture supplement, antigen of IVD reagents

**Image**

Bioactivity:



**For research use only**