

## Mouse PDGF-AA, Tag Free

---

Catalog	Unit Size
TBP0176-50	50 µg
TBP0176-100	100 µg

### Description

Platelet-derived growth factor subunit A is a human protein produced by the PDGFA gene. This protein belongs to the platelet-derived growth factor family.

Synonyms: Platelet-derived growth factor subunit A; PDGFAA; Platelet-Derived Growth Factor-AA; Glioma-derived growth factor; GDGF; Osteosarcoma-derived Growth Factor; ODGF

### Specifications

**Purity:** >95% as determined by SDS-PAGE.

**Endotoxin:** Less than 1.0 EU per µg by the LAL method.


### Formulation

Liquid in PBS (pH7.4), 5% trehalose, 5% mannitol

### Storage

Samples are stable for up to twelve months from date of receipt at -20°C to -80°C, avoid repeated freeze-thaw cycles.

### Data



Protein on SDS-PAGE under reducing (R) and non-reducing (NR) conditions.