

Oligo(dT)18 Primer

Catalog	Unit
TBS6120	500 μ l

Description

This product contains 500 μ l of 10 μ M Oligo(dT)18 primer, dissolved in ddH₂O.

Component

Product	Unit
Oligo(dT)18 Primer	500 μ l, 10 μ M

Specifications

Recommended Protocol: This product is intended for use in reverse transcription cDNA synthesis reactions, and should be added to reactions at a final concentration of 0.5-1 μ M.

Storage Condition: This product can be stored at -20°C for a maximum of 60 months. Avoid freeze thaw cycles.

For research use only