

Monoclonal Anti-Follicle Stimulating Hormone (FSH) Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10116-0.25MG	0.25mg
TBS10116-1MG	1mg

Description:	Mouse Monoclonal Anti-Follicle Stimulating Hormone (FSH) Antibody.
Specificity:	FSH 100% hCG 0.01% LH 0.01% TSH 0.02%
Host Animal:	Mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Human FSH.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	2.7 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10121 for human FSH detection.
Purity:	>95%.
Application:	For human FSH quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	-20°C for 2 years.

For Research Use only