
Specification Sheet

Catalog Number	Unit Size
TBS10160-0.2MG	0.2mg
TBS10160-1MG	1mg

Description:	Mab to methaqualone.
Specificity:	This antibody reacts with methaqualone.
Host Animal:	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Methaqualone-BSA.
Format:	Purified, liquid.
Conjugate:	No.
Purification:	DEAE chromatography.
Concentration:	3.9 mg/ml (OD280 nm, E ^{0.1%} = 1.40).
Affinity Constant:	Not determined.
Buffer:	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	N/A.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For methaqualone quantitative assays by EIA.
Storage:	2-8°C for 2 years. Long-term store at -20°C.

For Research Use only