

Catalog
TBS8054**Unit Size**
500 ml**Description**

Tribio Aspergillus Complete Growth Medium is a nutritionally rich, chemically defined, ready-to-use medium formulated for the cultivation and maintenance of *Aspergillus* species and other filamentous fungi. It supports rapid vegetative growth and efficient sporulation and is suitable for both routine culture and advanced research applications.

Application

- Routine cultivation of *Aspergillus* species
- Fungal genetics and mutant screening
- Sporulation studies
- Transformation recovery and strain maintenance
- General laboratory propagation of filamentous fungi

Main Features

- Ready-to-Use: No further preparation required.
- Supports rapid mycelial growth
- Promotes robust sporulation

Storage

Store at 4°C-8°C.

Relative Products

TBS2005	Protein Assay Kit
TBS5001	Protein Lysis Buffer
TBS5003-0	1xPBS
TBS5011-10	10 x PBST
TBS5008	TBS
TBS5008T	TBST
TBS5037	Urea Solution (8M)
TBS5038	Ammonium Sulfate Solution
TBS5040	0.5M EDTA Solution
TBS5041	MOPS Buffer
TBS5048	2% BSA in PBS
TBS5049-2	2% BSA in TBST
TBS5049-5	5% BSA in TBST
TBS5057	0.1% BSA ELISA Assay Buffer
TBS5077	Antigen retrieval Citrate Solution-pH6.0
TBS8024	MCF10 Cell Complete Medium
TBS8027	EMEM with L-Glutamine
TBS8030	McCoy's 5A Medium
TBS8031	30% BSA Solution in DPBS
TBS8047	MCF-7 Cell Complete Medium
TBS8056	LB Medium
TBS8057	SOB Medium
TBS8058	SOC Medium
TBS8059	2x YT Medium
TBS8060	Terrific Broth (TB) Medium
TBS8061	DMEM-high Glucose Medium

This product is for research use only.