

Monoclonal Anti-human Lactoferrin Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10146-0.2MG	0.2mg
TBS10146-1MG	1mg

Description:	Monoclonal anti-human lactoferrin antibody. This item (TBS10146, capture) can be paired with TBS10147 (detection) for ELISA.
Specificity:	This antibody reacts specifically with human lactoferrin (100%). It does not react with bovine lactoferrin.
Host Animal:	Balb/c mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Human lactoferrin.
Format:	Purified, liquid.
Conjugate:	No.
Purification:	DEAE chromatography.
Concentration:	2.6 mg/ml.
Affinity Constant:	Not determined.
Buffer:	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10147.
Purity:	>95%.
Application:	For human lactoferrin quantitative assays by EIA.
Storage:	2-8°C for 1 year. Long-term store at -20°C.

For Research Use only