

Monoclonal Anti-human Progesterone Antibody

Catalog Number	Unit Size
TBS10177-0.25MG	0.25mg
TBS10177-1MG	1mg

Description:	Mouse monoclonal anti-human progesterone antibody.
Specificity:	This antibody reacts specifically to human progesterone.
Host Animal:	Mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Progesterone-BSA.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	3 mg/ml.
Affinity Constant:	Not determined.
Buffer:	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH7.2
Pairing Assays:	N/A.
Purity:	>95%.
Application:	For human progesterone quantitative assays by EIA.
Storage:	Long-term store at -20°C.

For Research Use Only