

Catalog
TBS8128

Unit Size
500ml

Description

Tribio Human Diploid Cell Virus Maintenance Medium is a ready-to-use, liquid medium formulated to support the production of human diploid cells for viral applications, particularly rabies and hepatitis A viruses. It maintains the health and viability of human diploid cell lines during viral infection, providing a stable environment that promotes efficient virus replication.

Synonym

- HDC Virus Maintenance Medium.
- Human Diploid Cell Maintenance Medium (virus support).
- Diploid Cell Virus Support Medium (liquid).
- Virus-Compatible Diploid Maintenance Medium.

Application

- Maintenance of human diploid cells during virus infection/replication phases.
- Viral propagation workflows requiring stable cell viability post-infection.
- Routine virology studies in diploid cell systems.
- Bioprocessing applications where serum reduction or removal is required.

Properties

Concentration:	Ready-to-use
Form:	Liquid
Sterility:	Sterile filtered
Composition:	[+] L-glutamine, [+] phenol red, [-] sodium bicarbonate
Storage:	Store at 2-8°C.

Relative Products

TBS2005	Protein Assay Kit
TBS5001	Protein Lysis Buffer
TBS5003-0	1xPBS
TBS5011-10	10 x PBST
TBS5008	TBS
TBS5008T	TBST
TBS5037	Urea Solution (8M)
TBS5038	Ammonium Sulfate Solution
TBS5040	0.5M EDTA Solution
TBS5041	MOPS Buffer
TBS5048	2% BSA in PBS
TBS5049-2	2% BSA in TBST
TBS5049-5	5% BSA in TBST
TBS5057	0.1% BSA ELISA Assay Buffer
TBS5077	Antigen retrieval Citrate Solution-pH6.0
TBS8024	MCF10 Cell Complete Medium
TBS8027	EMEM with L-Glutamine
TBS8030	McCoy's 5A Medium
TBS8031	30% BSA Solution in DPBS
TBS8047	MCF-7 Cell Complete Medium
TBS8056	LB Medium
TBS8057	SOB Medium
TBS8058	SOC Medium
TBS8059	2x YT Medium
TBS8060	Terrific Broth (TB) Medium
TBS8061	DMEM-high Glucose Medium

This product is for research use only.