

Monoclonal Anti-17-Hydroxypregnenolone (17-OHP) Antibody

Catalog Number	Unit Size
TBS10406-0.2MG	0.2mg
TBS10406-1MG	1mg

Description:	Mouse monoclonal anti-17-Hydroxypregnenolone (17-OHP) antibody	
Specificity:	17-hydroxypregnenolone	100%
	17-hydroxyprogesterone	120%
	Pregnenolone	0.1%
	Progesterone	0.1%
Host Animal:	Mouse	
Isotype:	IgG2a	
Source:	Mouse ascites	
Immunogen:	17-hydroxypregnenolone-KLH	
Format:	Purified, liquid	
Purification:	DEAE chromatography	
Concentration:	4.59 mg/ml	
Affinity Constant:	Not determined	
Buffer:	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH7.2	
Pairing Assays:	N/A	
Purity:	>95%	
Application:	For 17-Hydroxypregnenolone quantitative assays by EIA	
Storage:	Long-term store at -20°C	

For Research Use Only