

**Catalog**  
TBS8133

**Unit Size**  
500ml

### Description

Tribio Serum-Reduced BHK-21 Medium is a specially formulated culture medium designed to support the **growth** and **maintenance** of BHK-21 (Baby Hamster Kidney) cells under low-serum conditions. Enriched with essential nutrients, vitamins, and salts optimized for BHK-21 cells, it promotes healthy proliferation while reducing the need for high serum concentrations. Supplied as a ready-to-use liquid, this medium delivers consistent performance and reliable results, making it ideal for serum-reduced **suspension** culture in both research and bioproduction settings.

### Application

- Routine culture of BHK-21 cells under low-serum conditions.
- Primarily designed for serum-reduced suspension cultivation of BHK-21 cells.
- Recombinant protein expression.
- Serum-reduction process optimization and scale-up.

### Properties

Concentration:	Ready-to-use
Form:	Liquid
Sterility:	Sterile filtered
Composition:	[+] L-glutamine, [+] Poloxamer 188, [-] Phenol red, [-] Sodium bicarbonate
Storage:	Store at 2 - 8°C

### Relative Products

TBS2005	<a href="#">Protein Assay Kit</a>
TBS5001	<a href="#">Protein Lysis Buffer</a>
TBS5003-0	<a href="#">1xPBS</a>
TBS5011-10	<a href="#">10 x PBST</a>
TBS5008	<a href="#">TBS</a>
TBS5008T	<a href="#">TBST</a>
TBS5037	<a href="#">Urea Solution (8M)</a>
TBS5038	<a href="#">Ammonium Sulfate Solution</a>
TBS5040	<a href="#">0.5M EDTA Solution</a>
TBS5041	<a href="#">MOPS Buffer</a>
TBS5048	<a href="#">2% BSA in PBS</a>
TBS5049-2	<a href="#">2% BSA in TBST</a>
TBS5049-5	<a href="#">5% BSA in TBST</a>
TBS5057	<a href="#">0.1% BSA ELISA Assay Buffer</a>
TBS5077	<a href="#">Antigen retrieval Citrate Solution-pH6.0</a>
TBS8024	<a href="#">MCF10 Cell Complete Medium</a>
TBS8027	<a href="#">EMEM with L-Glutamine</a>
TBS8030	<a href="#">McCoy's 5A Medium</a>
TBS8031	<a href="#">30% BSA Solution in DPBS</a>
TBS8047	<a href="#">MCF-7 Cell Complete Medium</a>
TBS8056	<a href="#">LB Medium</a>
TBS8057	<a href="#">SOB Medium</a>
TBS8058	<a href="#">SOC Medium</a>
TBS8059	<a href="#">2x YT Medium</a>
TBS8060	<a href="#">Terrific Broth (TB) Medium</a>
TBS8061	<a href="#">DMEM-high Glucose Medium</a>

**This product is for research use only.**