

Human BDNF Recombinant Protein (N-His)

Catalog	Unit
TBP0238-50	50 UG
TBP0238-100	100 UG
TBP0238-1	1 MG

Product Details

Description: Recombinant Human BDNF Protein, N-His expressed in E. coli, spanning His129-Arg247.

Expression system: E. coli

Accession: P23560

Protein length: His129-Arg247

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 17.52 kDa


Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: ProBDNF, Abrineurin, Brain-derived neurotrophic factor, BDNF

Image



SDS-PAGE for Human BDNF Recombinant Protein (N-His)

For research use only