

Mouse BDNF Recombinant Protein (N-His)

Catalog	Unit
TBP0239-50	50 UG
TBP0239-100	100 UG
TBP0239-1	1 MG

Product Details

Description: Recombinant BDNF Protein, N-His expressed in E. coli.

Expression system: E. coli

Accession: P23560, P21237, P23363, Q7YRB4, Q0EAB7, Q95106, P14082

Protein length: P23560[His129-Arg247], P21237[His131-Arg249], P23363[His131-Arg249], Q7YRB4[His129-Arg247], Q0EAB7[His129-Arg247], Q95106[His132-Arg250], P14082[His134-Arg252]

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Mouse, Rat, Dog, Horse, Bovine, Pig

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 15.96 kDa


Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: Brain-derived neurotrophic factor, BDNF, BDNF precursor form, ProBDNF, Bdnf

Image



SDS-PAGE for Recombinant Mouse BDNF Protein, N-His

For research use only