

Human FGF9 Recombinant Protein (N-Fc)

Catalog	Unit
TBP0241-50	50 UG
TBP0241-100	100 UG
TBP0241-1	1 MG

Product Details

Description: Recombinant Human FGF9 Protein, N-Fc expressed in Mammalian Cells, spanning Leu4-Ser208.

Expression system: Mammalian Cells

Accession: P31371

Protein length: Leu4-Ser208

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >95%

Predicted molecular weight: 50.97 kDa


Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: FGF9, Glia-activating factor, GAF, Heparin-binding growth factor 9, FGF-9, Fibroblast growth factor 9, HBGF-9

Image



SDS-PAGE for Recombinant Human FGF9 Protein, N-Fc

For research use only