

Human SEMA3A Recombinant Protein (N-His)

Catalog	Unit
TBP0243-50	50 UG
TBP0243-100	100 UG
TBP0243-1	1 MG

Product Details

Description: Recombinant Human SEMA3A Protein, N-His expressed in E. coli, spanning Lys26-His276.

Expression system: E. coli

Accession: Q14563

Protein length: Lys26-His276

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 31.42 kDa

Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: Semaphorin III, Sema III, SEMA3A, Semaphorin-3A, SEMAD

For research use only