

Human NRG2 Recombinant Protein (N-His)

Catalog	Unit
TBP0244-50	50 UG
TBP0244-100	100 UG
TBP0244-1	1 MG

Product Details

Description: Recombinant Human NRG2 Protein, N-His expressed in E. coli, spanning Ser114-Asn374.

Expression system: E. coli

Accession: O14511

Protein length: Ser114-Asn374

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 31.02 kDa

Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: Pro-neuregulin-2, membrane-bound isoform, Pro-NRG2, Neuregulin-2, NRG-2, Divergent of neuregulin-1, DON-1, Neural- and thymus-derived activator for ERBB kinases, NTAK, NRG2, NTAK

Image



SDS-PAGE for Human NRG2 Recombinant Protein (N-His)

For research use only