

Human SLIT3 Recombinant Protein (N-His)

Catalog	Unit
TBP0245-50	50 UG
TBP0245-100	100 UG
TBP0245-1	1 MG

Product Details

Description: Recombinant Human SLIT3 Protein, N-His expressed in E. coli, spanning Ala33-Val262.

Purity: >90%

Concentration: 0.56 mg/mL

Expression system: E. coli

Purification: Affinity-chromatography

Accession: O75094

Protein length: Ala33-Val262

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant


Predicted molecular weight: 28.68 kDa

Formulation: Liquid in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: Slit homolog 3 protein, Multiple epidermal growth factor-like domains protein 5, MEGF5, SLIL2, Slit-3, SLIT3, Multiple EGF-like domains protein 5, KIAA0814

Image



SDS-PAGE for Human SLIT3 Recombinant Protein (N-His)

For research use only