

Mouse SLIT3 Recombinant Protein (N-His)

Catalog	Unit
TBP0246-50	50 UG
TBP0246-100	100 UG
TBP0246-1	1 MG

Product Details

Description: Recombinant Mouse SLIT3 Protein, N-His expressed in E. coli, spanning Ala33-Val262.

Expression system: E. coli

Accession: Q9WVB4

Protein length: Ala33-Val262

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Mus musculus (Mouse)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 28.68 kDa

Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: Slit homolog 3 protein, Multiple epidermal growth factor-like domains protein 5, MEGF5, SLIL2, Slit-3, SLIT3, Multiple EGF-like domains protein 5, KIAA0814

For research use only