

Human NRG3 Recombinant Protein (N-GST & C-His)

Catalog	Unit
TBP0247-50	50 UG
TBP0247-100	100 UG
TBP0247-1	1 MG

Product Details

Expression system: E. coli

Accession: P56975

Protein length: Glu282-Pro342

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 35.2 kDa

Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: NRG-3, NRG3, Neuregulin-3, Pro-NRG3, Pro-neuregulin-3, membrane-bound isoform

For research use only