

Human FLRT3 Recombinant Protein (N-His)

Catalog	Unit
TBP0248-50	50 UG
TBP0248-100	100 UG
TBP0248-1	1 MG

Product Details

Expression system: E. coli

Accession: Q9NZU0

Protein length: Ser30-Asn349

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 38.82 kDa

Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: FLRT3, Fibronectin-like domain-containing leucine-rich transmembrane protein 3, KIAA1469, Leucine-rich repeat transmembrane protein FLRT3

For research use only