

Monoclonal Anti-human D-dimer Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10109-0.2MG	0.2mg
TBS10109-1MG	1mg

Description:	Monoclonal Anti-human D-dimer antibody. This item (TBS10109) can be paired with TBS10110 for ELISA.
Specificity:	This antibody reacts specifically with human D-Dimer (100%), does not react with human fibrinogen.
Host Animal:	Balb/c mouse.
Source:	Mouse ascites.
Immunogen:	Human D-Dimer.
Format:	Purified, liquid.
Purification:	Protein G chromatography.
Concentration:	3.5 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10110 for human D-Dimer detection.
Purity:	>95%.
Application:	For human D-Dimer quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	2-8°C for 1 year. Long-term store at -20°C.

For Research Use only