

Anti-Human Ganglioside GD2 Mouse Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10258-0.1	0.1 MG
TBS10258-0.5	0.5 MG
TBS10258-1	1 MG

Description:	Monoclonal Anti-human Ganglioside GD2 antibody.
Species reactivity:	Human
Concentration:	3 mg/mL
Form:	Liquid
Storage buffer:	0.01M PBS pH 7.4
Purity:	>95%
Clonality:	Monoclonal
Applications:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo.
Target:	Ganglioside GD2
Purification:	Protein A/G purified from cell culture supernatant.
Stability and Storage:	Store at -20 to -80°C for twelve months from the date of receipt.

For Research Use only