

Anti- Human DCX Rabbit Polyclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10519-100	100 UG

Description:	Neuronal migration protein doublecortin (DCX) is a microtubule-associated protein with a molecular weight of approximately 40 kDa. It plays an essential role in the early stages of neuronal migration and cortical layer formation during cerebral cortex development.
Host species:	Rabbit
Species reactivity:	Human
Immunogen:	E. coli - derived recombinant Human DCX
Form:	Liquid
Storage buffer:	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Clonality:	Polyclonal
Applications:	ELISA, IHC, WB
Target:	Lis-X, LISX, DCX, Neuronal migration protein doublecortin, Dublin, Lissencephalin-X, DBCN
Purification:	Purified by antigen affinity column.
Stability and Storage:	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

For Research Use only