

Catalog	Unit
TBP0169-0.1	0.1 mg
TBP0169-1	1 mg
TBP0169-5	5 mg
TBP0169-10	10 mg

Description

Tribioscience Recombinant Human Insulin-like Growth Factor (rhIGF-1) is manufactured through a recombinant expression system, guaranteeing consistent quality and high purity. This animal-free protein, with its low endotoxin levels, offers superior safety and reliability for research use.

Specification

Source: Escherichia coli

Appearance: White lyophilized powder

Molecular weight: 8.3 kDa

Formulation: It is lyophilized with mannitol as the stabilizer

Purity: More than 95% as determined by SDS-PAGE

Endotoxin: Less than 5 EU/mg


Bioactivity: $\geq 1.0 \times 10^6$ IU/mg

Storage: Stored at temperatures -20°C for at least 36 months

Applications: Cell culture supplement, antigen of IVD reagents

Image

Bioactivity:



For research use only