

ELISA Dilution Buffer (0.1% BSA in TBST) (TBS5057)

Catalog
TBS5057-500

Unit Size
500 ml

DESCRIPTION

Tribio ELISA Dilution Buffer is a ready-to-use buffer containing 0.1% BSA in TBST, designed for dilution of ELISA samples and antibodies. It helps support protein stability and reduce nonspecific binding and background.

APPLICATION

- ELISA sample dilution
- Antibody dilution
- ELISA assay procedures
- Protein and antibody handling

PROPERTIES

Concentration:	Ready-to-use
Form:	Liquid
Sterility:	Sterile filtered
Composition:	0.1% BSA in TBST
Storage:	Store at 2-8°C

SHELF LIFE

12 months from the date of manufacture when stored under recommended conditions.

PROCEDURE

Mix gently before use. Use directly as supplied according to the specific ELISA assay protocol. Avoid contamination during handling.

RELATIVE PRODUCTS

TBS2005	Protein Assay Kit
TBS5001	Protein Lysis Buffer
TBS5003-0	1xPBS
TBS5011-10	10 x PBST
TBS5008	TBS
TBS5008T	TBST
TBS5037	Urea Solution (8M)
TBS5038	Ammonium Sulfate Solution
TBS5040	0.5M EDTA Solution
TBS5041	MOPS Buffer
TBS5048	2% BSA in PBS
TBS5049-2	2% BSA in TBST
TBS5049-5	5% BSA in TBST
TBS5057	0.1% BSA ELISA Assay Buffer
TBS5077	Antigen retrieval Citrate Solution-pH6.0
TBS8024	MCF10 Cell Complete Medium
TBS8027	EMEM with L-Glutamine
TBS8030	McCoy's 5A Medium
TBS8031	30% BSA Solution in DPBS
TBS8047	MCF-7 Cell Complete Medium
TBS8056	LB Medium
TBS8057	SOB Medium
TBS8058	SOC Medium
TBS8059	2x YT Medium
TBS8060	Terrific Broth (TB) Medium
TBS8061	DMEM-high Glucose Medium

This product is for research use only.