

Human ETNK1 Recombinant Protein (N-His) Bioactive

Catalog	Unit
TBP0249-50	50 UG
TBP0249-100	100 UG
TBP0249-1	1 MG

Product Details

Expression system: E. coli

Accession: Q9HBU6

Protein length: Ile94-Glu452

Applications: ETNK1 enzyme Activity, ELISA, Immunogen, SDS-PAGE, WB. ETNK1 bioactivity was validated by ETNK enzyme activity, as shown in Fig. 1.

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 43.77 kDa

Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: EC:2.7.1.82, EKI 1, EKI1, ETNK1, Ethanolamine kinase 1

For research use only

