

Human CS Recombinant Protein (N-His)

Catalog	Unit
TBP0253-50	50 UG
TBP0253-100	100 UG
TBP0253-1	1 MG

Product Details

Expression system: E. coli

Accession: O75390

Protein length: Lys103-Ala394

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 34.9 kDa


Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: CS, Citrate synthase, mitochondrial, Citrate (Si)-synthase

SDS-PAGE Analysis Result



For research use only