

Human RB1 Recombinant Protein (N-GST)

Catalog	Unit
TBP0254-50	50 UG
TBP0254-100	100 UG
TBP0254-1	1 MG

Product Details

Expression system: E. coli

Accession: P06400

Protein length: Gln639-Thr778

Applications: ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress

Species: Homo sapiens (Human)

Nature: Recombinant

Purity: >90%

Predicted molecular weight: 43.67 kDa


Form: Lyophilized

Storage buffer: Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1 mM EDTA, 4% Trehalose, 1% Mannitol.

Stability and Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Alternate Names: pRb, p110-RB1, Retinoblastoma-associated protein, Rb, p105-Rb, RB1, pp110

SDS-PAGE Analysis Result



For research use only