

Specification Sheet

Catalog Number	Unit Size
TBS10002-0.2MG	0.2mg
TBS10002-1MG	1mg

Description:	Mab to benzodiazepine.
Specificity:	This antibody reacts with benzodiazepine.
Host Animal:	Mouse.
Isotype:	IgG1, kappa light chain.
Source:	Mouse ascites.
Immunogen:	Benzodiazepine-KLH.
Format:	Purified, liquid.
Conjugate:	No.
Purification:	DEAE chromatography.
Concentration:	3.60 mg/ml.
Affinity Constant:	Not determined.
Buffer:	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Not determined.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For benzodiazepine quantitative assays by EIA.
Storage:	2-8°C for 2 years. Long-term store at -20°C.

For Research Use only