

Anti-Mullerian Hormone (Anti-AMH) Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10013-0.25MG	0.25mg
TBS10013-1MG	1mg

Description:	Mouse Monoclonal Antibody to Anti-mullerian Hormone (Anti-AMH)
Specificity:	This antibody is specific to human AMH
Host Animal:	Mouse
Isotype:	IgG1
Source:	Mouse ascites
Immunogen:	Human AMH, recombinant from HEK 293 cells
Format:	Purified, liquid
Purification:	DEAE chromatography.
Concentration:	4.20 mg/ml
Affinity Constant:	Not determined
Buffer:	0.015M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH 7.2
Pairing Assays:	Pairs with and TBS10012
Purity:	>95% by HPLC and SDS-PAGE
Application:	For human AMH quantitative assays by EIA
Storage:	2-8°C for 1 year. Long-term store at -20°C

For Research Use only